



ENDOSCOPE CLEANING GUIDE

EVIS™ 40/140/240 & EXERA™ 160-Series GI Endoscopes

Refer to companion *ENDOSCOPE DISINFECTION GUIDE* for disinfection information

WARNING: This guide is only a summary of the steps necessary to properly clean your endoscope. Be sure to follow the detailed steps outlined in the *ENDOSCOPE REPROCESSING MANUAL* that was included with your endoscope when purchased.

PRECLEANING - IMMEDIATELY AFTER USE

1

Wipe down the insertion tube with a detergent-soaked cloth.

2

Aspirate detergent, then air through the suction channel.

3

Install AW channel-cleaning adapter. Flush water, then air, through air/water channels.

4 For endoscopes with elevators or auxiliary-water feeding only:

Flush detergent, then air, into elevator-wire channel/auxiliary-water channel.

5

Disconnect all detachable parts and cleaning equipment.

LEAKAGE TESTING

6

Attach water-resistant cap (videoscopes only).

7 **CAUTION:** If a leak is detected, the endoscope must be repaired. Follow manufacturer's instructions.

Attach leakage tester. Turn on pump. Prior to immersion, confirm bending section has expanded. Immerse entire endoscope in clean water. Perform leakage test. Angulate tip during test. Remove from water. Detach tester after tip has deflated.

MANUAL CLEANING

8

Immerse in freshly prepared detergent solution. Clean all external surfaces. Brush endoscope distal tip (and forceps-elevator, if applicable).

9

Brush the insertion tube portion of the suction channel. Repeat until all debris is removed.

10

Brush the universal cord portion of the suction channel. Repeat until all debris is removed.

11

Brush the suction cylinder and instrument channel port. Repeat until all debris is removed.

12

Aspirate detergent solution into suction channel.

13

Attach injection tube and channel plug. Flush detergent solution into air, water and suction channels. Brush forceps-elevator, if applicable.

14 For endoscopes with elevators or auxiliary-water feeding only:

Flush detergent solution into auxiliary-water/elevator-wire channel.

15

Soak in detergent solution for the time recommended by the detergent manufacturer.

RINSING

16

Rinse the entire instrument in clean water.

17

Flush clean water through all channels (including auxiliary-water/elevator-wire, if applicable). Flush air through all channels.

18

Use a soft, lint-free cloth to remove excess moisture from endoscope and cleaning accessories in preparation for disinfection.

IMPORTANT

- Wear all appropriate personal protective equipment.
- Meticulous cleaning is essential for effective disinfection.
- Be sure to reprocess all removable parts (e.g., valves) according to the reprocessing manual.

Follow cleaning with appropriate high-level disinfection and alcohol-air drying.