

INSTRUCTION FOR USE

Instrument Cleaner; Concentrated



PRODUCT NAME : Instrument Cleaner; Concentrated

TRADEMARK : LIQUISEPT™ QAS CON

PRODUCT REF NO: 6018

PACKING TYPE	ORDER NO
500 ml HDPE; Dosing Bottle	436018140500
1 Liter HDPE; Dosing Bottle	436018131000
5 Liter HDPE; Canister	436018025000

INTENDED USE: Used for cleaning and pre-disinfection of medico-surgical instruments, medical devices, thermosensitive instruments endoscopes. Suitable for soaking baths, ultrasonic bins and automatic washing machines.

PRODUCT SPECIFICATION

Form	Liquid
Color	Yellowish
Odor	Specific
pH @ 25°C (ca):	7.5±0.3

COMPOSITION

Didecyl dimethyl ammonium chloride (6%), alkyl dimethyl benzyl ammonium chloride (5%), isopropanol (5%), detergents, chelating agents, defoamer, thickening agents, emulsion agents, protease, buffers, corrosion inhibitors and stabilizing agent. Without aldehyde and phenol.

MICROBIOLOGICAL PROPERTIES

Effect	Standard	Microorganism (with ATCC)	Concentration of use	Contact Time
Bactericidal	EN 13727 / EN 14561	<i>Staphylococcus aureus</i> ATCC 6538	1/100	15 min.
		<i>Pseudomonas aeruginosa</i> ATCC 15442		
		<i>Enterococcus hirae</i> ATCC 10541		
Fungicidal	EN 13624 / EN 14562	<i>Candida albicans</i> ATCC 10231	1/100	15 min.
		<i>Aspergillus brasiliensis</i> ATCC 16404		60 min.

*The product is the activity table with the specified concentration and duration of use.

HOW TO USE

Dilute concentrated product with preferably distilled or tap water with a 1/100 ratio. It means for nearly 2 liters (1.980,00 ml) of water, you must add 20 ml (1 pump) of a concentrated product. Immerse instruments inside a soaking

bath full of ready-to-use products for 15 minutes. Rinse preferably with sterile distilled water. Prepare new solutions every day.

WARNINGS AND POTENTIAL ADVERSE EFFECTS

Causes severe skin burns and eye damage. Do not breathe dust/fumes/gas/mist/vapors/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF INHALED: Remove the person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do continue rinsing.

STABILITY AND STORAGE

Product; in the shade, in the original packaging, the package should be kept tightly closed and in a well-ventilated environment.

Storage Temperature	Shelf Life	Usage Time After Opening the Cap
5/25 °C	2 years	3 months

PACKAGING

Product	Packaging Type	Packaging Units
LIQUISEPT™ QAS CON	500 ml HDPE; Dosing Bottle	500 ml x 24
LIQUISEPT™ QAS CON	1 Liter HDPE; Dosing Bottle	1L x 12
LIQUISEPT™ QAS CON	5 Liter HDPE; Canister	5L x 2

LABEL SYMBOL KEYS



Catalog
Number



Production
Date



Storage
Temperature



Manufacturer
Information



Batch Code



Expiry Date



Refer to User
Guide



Recycling
Information



GHS05
Corrosive

MANUFACTURER

GBL® Gül Biyoloji Laboratuvarı Sanayi ve Ticaret Anonim Şirketi

HQ: Serifali Mah. Hattat Sk. No:10 P.O.: 34775 Ümraniye İstanbul TÜRKİYE

Factory: Dudullu OSB Mah. İMES C Blok 305 Sk. No:16 P.O.: 34775 Ümraniye İstanbul TÜRKİYE

Phone: +90 216 364 15 00 **Email:** export@gbl.com.tr **Web:** www.gbl.com.tr