



6106 - CITRO PLUS

Hemodialysis Machine Disinfectant; Citric Acid 21%



INDICATION

Cleaning, thermochemical disinfection, decalcification, removal of inorganic and organic residues in hemodialysis machines.

COMPOSITION

Citric acid (21%), Lactic acid (2.5%), Malic acid (2.5%), buffers, preservatives and RO water.

MICROBIOLOGICAL PROPERTIES

Bactericide, fungicide, mycobactericide, virucide, sporicide.

PACKAGING








1 carton: 5 liter can x 2

1 pallet: cartons / cans

Order .no: 436106035000

INSTRUCTION

CITRO PLUS™ is a potent disinfectant solution because of its components synergistic effect. CITRO PLUS™ at 84°C has a broad spectrum of microbicidal activity. CITRO PLUS™ at 37°C decalcifies very effectively organic and inorganic residues. Product to be diluted 1/25 automatically by the machine. Please refer to the hemodialysis machine's "operating manual" for more details.

-  Dissolves blood residue.
-  Perfectly removes CaCO₃.
-  Disinfection and decalcification in one process.
-  Biodegradable, non-toxic.
-  Active ingredients composed of natural substances.
-  Colorless; almost odorless.
-  Safe for health and environment.